

## Technical Information Piece

Date: April 11, 2016  
 To: All Field Personnel  
 From: BIOTRONIK Technical Services  
 Re: Indication Specific Program Suggestions for BioMonitor 2

---

BioMonitor 2 allows the user to monitor the patient for various clinical events. The BioMonitor 2 has five automatic triggers and a patient trigger option to capture these events. To help guide programming, the table below provides programming suggestions for different indications. This is designed to make programming easier and tailor the device to meet the clinical needs of the patient while reducing non-actionable recordings.

Parameter	Default Setting	Syncope	Palpitations	AF Management	Suspected AF	Post AF Ablation
AF						
AF Sensitivity	Medium	Low	Medium	Medium	Low	Medium
Confirmation Time (minutes)	6	10	6	2	6	2
Bigeminy Rejection	Standard	Aggressive	Standard	Standard	Standard	Standard
High Ventricular Rate (HVR)						
HVR Limit	180	180	180	180	180	180
HVR Counter	16	32	32	32	32	32
Bradycardia						
Brady Zone Limit	40	30	30	30	30	30
Brady Duration	10	10	30	30	30	30
Asystole Duration	3	3	3	OFF	OFF	OFF
Sudden Rate Drop	OFF	On	OFF	OFF	OFF	OFF

Decrease (%)		50				
Sensitivity		Low				

To save time programming devices, the programmer can store up to three program sets under the parameter page of the device. To set up predetermined programs, do the following:

1. Program the desired settings

Diagnostics		Home Monitoring		Patient	
<b>Detection settings</b>			<b>Sensing settings</b>		
Atrial fibrillation (AF)	<b>Individual</b>	SensingConsult	<b>PVC detection</b>		
High ventricular rate (HVR)	<b>ON</b>				
Bradycardia	<b>ON</b>				
Sudden rate drop (SRD)	<b>ON</b>				
Asystole duration [s]	<b>3</b>				
Patient trigger	<b>ON</b>				
<b>Resting rate period</b>					
Start resting period [hh:mm]	<b>02:00</b>				
Resting period duration [h]	<b>4.0</b>				
		<b>Help</b>	<b>Program sets</b>	<b>Program</b>	

2. Press the button **Program sets**
3. Press the **Save** button. A keyboard will appear allowing the user to type in the name of the program. Once the name is typed in, press the enter key.
4. Repeat until three programs are entered

### Program sets

#### Show standard program

Show first interrogated program

Show permanent program

#### Syncopy

Save

#### AF Management

Save

#### Suspected AF

Save

When programming a device, the user can select one of the program sets to populate the parameters. A message, **Stored program displayed**, will appear when one of the program sets is selected.

Stored program displayed

PSA

Diagnostics

Home Monitoring

Patient

Detection settings

Atrial fibrillation (AF)

Individual

High ventricular rate (HVR)

ON

Bradycardia

ON

Sudden rate drop (SRD)

ON

Asystole duration [s]

3

Patient trigger

ON

Sensing settings

SensingConsult

PVC detection

Resting rate period

Start resting period [hh:mm]

02:00

Resting period duration [h]

4.0

Help

Program sets

Program

Press the **Program** button to transmit the program to the device.

Please contact Technical Services with any questions: 1-800-284-6689